



W-SMARTWAVE

High Energy Shockwave Therapy System

W-Smartwave Shockwave therapy system is a multidisciplinary device used in orthopedics, physiotherapy, sports medicine, urology and veterinary medicine. Together with being a non-surgical therapy with no need for painkillers makes it an ideal therapy to speed up recovery and cure various indications causing acute or chronic pain

W-Smartwave Shockwave therapy is a method for the treatment of superficial orthopedic disorders in which the hand piece generates an electromagnetic pressure wave, which is transferred to the human body by an applicator head and which radiates out in the tissue.

The electromagnetic energy is absorbed by the tissue and decreases as the distance from the application site increases. The mechanical stress triggers reactions in the tissue that have a positive influence on a wide range of orthopedic and neurological conditions. W-Smartwave is used for the following conditions:

- Insertion tendinopathy
- Tendon problems
- Chronic inflammation
- Hypertonic muscle systems
- Myofascial syndromes
- Hemotoma treatments



SPECIFICATIONS:

Display : Coloured 7 inch Touch Screen

Frequency : 1 to 22 Hz in steps of 0,5Hz

Energy : 10 to 190mJ (0-6bar)

Applicator life : 20 million cycles

Stimulation counter : Online and also storage facility

Types of Titanium Heads : 1238cm², 598cm², 278cm², 210cm²



Applicator: Latest Electromagnetic Technology with Forced air cooling Fan

Remote operation: Foot Switch